

How to Use An EpiPen®:

(**Note: up to 90% of the epinephrine liquid stays in the auto-injector once it has been used)




1. Pull off gray safety cap to activate




2. Hold EpiPen® in the middle – do not touch either end with your hand

3. Place the black tip firmly against the outer middle part of the thigh only and push in until you feel it click – remember: you can use EpiPen® over clothing in an emergency



4. Hold in place for at least 10-15 seconds then remove it



5. Call 911 immediately

6. Rub the site of the injection gently for at least 10 seconds.